




Age Prevention SX-4500

With young, fresh, healthy skin Neurotris advanced micro-current treatments with dynamic intelligence can slow the aging process, help maintain a youthful glow, deliver hydration balance and support cellular metabolism. The following treatments are based on skin types. If your client has minimal dysfunction like clogged pores or a few blemishes these treatment steps will apply.

- 🌿 Pre-Cleanse to remove cosmetics if applicable.
- 🌿 2nd Cleanse using water base foaming cleanser and utilize your Photon Probe to deep pore cleanse in the following manner.
 - Apply Positive Clip Cable to upper arm, insert photon probe. Choose Oxygen or Refine Nano and set amplitude to client comfort, set ultra-sonic on high for normal to dry skin and low for combination to oily skin. Set LED to Green for blemishes and Blue for oily to combination Skin. Use circular motions for 2-3 minutes covering the entire face, focus on blackheads if needed. Remove cleanser.
- 🌿 Gently exfoliate with enzymes, soft scrubs, or derma plane.
- 🌿 Perform extractions if needed no more than 3-5 minutes.
- 🌿 Using Infinity Probes large on black and small on red, apply Neurotris Conducting Gel,
 - Set to Refine Pico and adjust amplitude to client comfort. With probes set to 1 to 2 inches apart begin on the right and gently glide the probes across the neck and up to the chin, continue to glide in scissor like fashion over mouth, around the cheeks. Stopping under the eye. Station the negative probe then use the positive probe in circular motions in the corners of the

nose, glide over the cheeks continue gliding around the eyes and over the forehead and down the left side of the face in the same way you did the right. Eyes, cheek, nose, mouth, chin and neck. This total gentle gliding should take 3-5 minutes. Take care not to touch the tips of the Infinity probes and remember to keep the black above the red spaced no more than an inch apart.

 Use Neurotris Hydrating mist to rehydrate Neurotris conducting gel.

- Set the machine to Refine Pico adjusting amplitude to client comfort. Using the Infinity Probes large probe on black and small probe on red, focus on lips. Begin at the lower lip by setting the negative probe on lower lip fullness, pinching the lip between the probes concentrate on the perimeter of the lip holding 5-8 seconds as you go around the mouth for 2-3 minutes.
- Then move to the nasal lip folds with negative above positive lift and hold the mouth lines for 8-10 seconds do this on both sides of the mouth for 2-3 minutes.
- Keeping the negative about positive use the erase movement slowly over upper lip and on lip folds around the entire mouth. Do this for 1-2 minutes.
- Continue to the eye area holding the negative probe on the upper orbital bone, with circular motions apply the positive probe on the lower orbital bone and erase from the outside of the eye inwards for 1-3 minutes on each eye and in between the eyes on the glabella folds. Adjust amplitude to client comfort.

- 🌀 **Option A:** Refresh Neurotris conducting gel and switch to Gloves and clip cables.
 - Set to Tone begin using firm pressure and glide, lift and hold for 8-10 seconds at a time for each area of the neck, chin, cheeks, and forehead. Cover as much of the skin as possible taking care not to touch the gloves together. Do this for a total of 10-15 minutes, gently feather off.
- 🌀 **Option B:** Refresh Neurotris conducting gel and change to Pico Probes
 - Set to Tone using firm pressure glide, lift and hold for 8-10 seconds at a time on the neck, chin, cheeks, and forehead. Cover as much of the skin as possible taking care not to touch the Pico Probes together. Remember to keep the negative Pico Probe above the positive Pico probe, begin on the right and move down the left mirroring your movements. Do this for 10 -15 minutes in total.
- 🌀 Remove any traces of Neurotris gel and choose resting mask based on skin type. (cream for normal to dry, gel for combination or clay for oily or in T zones) remove as directed and tone if desired.
- 🌀 Apply water-based serums allow to absorb, apply a small amount of Neurotris conductive gel and penetrate product further by using Photon Probe (Nano Channel 1 Plugin) with the positive Infinity probe and small probe wand.
 - Set Ultrasonic to low, LED to green on the Photon probe, Set to Ion EMTA, adjust amplitude and use slow circular motions for 2-3 minutes covering the face and neck. Remove any conducting gel residue.
- 🌀 Finish with any additional serum, eye care, lip care, spf for day or nutrition/treat creams by night.
- 🌀 Final tip remember to set each channel at the same amplitude or intensity and document your settings.